

# SECTM1 antibody



Catalog Number: 107523

## Product name

SECTM1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SECTM1 Mouse Monoclonal antibody. Positive FC detected in K-562 cells. Positive WB detected in PC-3 cells, K-562 cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 27 kDa

## Preparation

This antibody was obtained by immunization of SECTM1 recombinant protein (Accession Number: XM\_005256392). Purification method: Caprylic acid/ammonium sulfate precipitation.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgM

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107523(SECTM1 Antibody) at dilution of 1:1000



1X10<sup>6</sup> K-562 cells were stained with 0.2ug SECTM1 antibody (Catalog No:107523, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1000.