

Zyxin antibody

Catalog Number: 107549

Product name

Zyxin antibody

Immunogen

[Human ZYX Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

Zyxin Mouse Monoclonal antibody. Positive WB detected in HeLa cells, A431 cells, MCF-7 cells. Positive IF detected in HeLa cells. Positive IHC detected in human kidney tissue, human breast cancer tissue, human cervical cancer tissue, human placenta tissue, human testis tissue. Observed molecular weight by Western-blot: 78kd

Preparation

This antibody was obtained by immunization of Zyxin recombinant protein (Accession Number: NM_001010972). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107549(ZYX antibody) at dilution of 1:1000




Immunohistochemistry of paraffin-embedded human kidney slide using Catalog No:107549(ZYX Antibody) at dilution of 1:50

Zyxin antibody

Catalog Number: 107549



Immunohistochemistry of paraffin-embedded human kidney slide using Catalog No:107549(ZYX Antibody) at dilution of 1:50



Immunofluorescent analysis of HeLa cells using Catalog No:107549(ZYX Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Mouse IgG