

SMARCD1 antibody



Catalog Number: 107589

Product name

Monoclonal

SMARCD1 antibody

Ig Type

Specificity

Mouse IgG1

Human; other species not tested.

Applications

Antibody description

ELISA, WB

SMARCD1 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, HepG2 cells, human brain tissue. Observed molecular weight by Western-blot: 58 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of SMARCD1 recombinant protein (Accession Number: BC009368). Purification method: Protein G purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107589(SMARCD1 antibody) at dilution of 1:1000