

# SMARCD1 antibody



Catalog Number: 107589

## Product name

SMARCD1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SMARCD1 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, HepG2 cells, human brain tissue. Observed molecular weight by Western-blot: 58 kDa

## Preparation

This antibody was obtained by immunization of SMARCD1 recombinant protein (Accession Number: BC009368). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107589(SMARCD1 antibody) at dilution of 1:1000