

# SMN (Human,Mouse,Rat) antibody



Catalog Number: 107593

## Product name

SMN (Human,Mouse,Rat) antibody

## Specificity

Human,Mouse,Rat; other species not tested.

## Antibody description

SMN (Human,Mouse,Rat) Mouse Monoclonal antibody. Positive WB detected in 293 cells, HepG2 cells, Raji cells. Observed molecular weight by Western-blot: 38 kDa

## Preparation

This antibody was obtained by immunization of SMN (Human,Mouse,Rat) recombinant protein (Accession Number: NM\_022875). Purification method: Protein A purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



WB results of Catalog No:107593(anti-SMN2).



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107593(SMN2 antibody) at dilution of 1:1000