

# SMN-Exon7 antibody



Catalog Number: 107595

## Product name

SMN-Exon7 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SMN-Exon7 Mouse Monoclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, HepG2 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of SMN-Exon7 recombinant protein (Accession Number: NM\_000344). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107595(SMN1 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:107595(SMN1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)