

# SNRPD2 antibody



Catalog Number: 107599

---

## Product name

SNRPD2 antibody

## Specificity

Human; other species not tested.

## Antibody description

SNRPD2 Mouse Monoclonal antibody. Positive WB detected in Hela cells, MCF7 cells. Observed molecular weight by Western-blot: 16kDa

## Preparation

This antibody was obtained by immunization of SNRPD2 recombinant protein (Accession Number: NM\_001384647). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107599(SNRPD2 antibody) at dilution of 1:1000