

ST6GALNAC6 antibody



Catalog Number: 107604

Product name

ST6GALNAC6 antibody

Specificity

Human; other species not tested.

Antibody description

ST6GALNAC6 Mouse Monoclonal antibody.
Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 38-42 kDa

Preparation

This antibody was obtained by immunization of ST6GALNAC6 recombinant protein (Accession Number: XM_011518610). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:107604(ST6GALNAC6 Antibody) at dilution of 1:500