

# ST6GALNAC6 antibody



Catalog Number: 107604

## Product name

Monoclonal

ST6GALNAC6 antibody

## Ig Type

## Specificity

Mouse IgG1

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

ST6GALNAC6 Mouse Monoclonal antibody.  
Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 38-42 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of ST6GALNAC6 recombinant protein (Accession Number: XM\_011518610). Purification method: Protein G purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:107604(ST6GALNAC6 Antibody) at dilution of 1:500

## Clonality