

STMN3 antibody

Catalog Number: 107609

Product name

STMN3 antibody

Specificity

Human; other species not tested.

Antibody description

STMN3 Mouse Monoclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of STMN3 recombinant protein (Accession Number: NM_015894). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107609(STMN3 antibody) at dilution of 1:1000