

Syntaxin 6 antibody

Catalog Number: 107615

Product name

Syntaxin 6 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Syntaxin 6 Mouse Monoclonal antibody. Positive IHC detected in human brain tissue. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 29kd

Preparation

This antibody was obtained by immunization of Syntaxin 6 recombinant protein (Accession Number: NM_005819). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107615(STX6 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human brain using Catalog No:107615(STX6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:107615(STX6 antibody) at dilution of 1:50 (under 40x lens)