

TMEM70 antibody

Catalog Number: 107628

Product name


TMEM70 antibody

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Validations




Antibody description

TMEM70 Mouse Monoclonal antibody. Positive IHC detected in human liver tissue, human pancreas tissue. Positive WB detected in A2780 cells. Observed molecular weight by Western-blot: 18 kDa

A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:107628(TMEM70 antibody) at dilution of 1:500

Preparation

This antibody was obtained by immunization of TMEM70 recombinant protein (Accession Number: NM_017866). Purification method: Protein A purified.




Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Immunohistochemical of paraffin-embedded human liver using Catalog No:107628(TMEM70 antibody) at dilution of 1:50 (under 40x lens)

Storage

Store at -20°C. DO NOT ALIQUOT



Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human liver using Catalog No:107628(TMEM70 antibody) at dilution of 1:50 (under 40x lens)