

# TP63 antibody



Catalog Number: 107632

---

## Product name

TP63 antibody

## Specificity

Human; other species not tested.

## Antibody description

TP63 Mouse Monoclonal antibody. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of TP63 recombinant protein (Accession Number: NM\_003722). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107632(TP63 Antibody) at dilution of 1:1000