

TRAF3 antibody



Catalog Number: 107636

Product name

TRAF3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TRAF3 Mouse Monoclonal antibody. Positive WB detected in Raji cells, LPS treated U-937 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of TRAF3 recombinant protein (Accession Number: XM_017021618). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:50-1:500

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107636(TRAF3 Antibody) at dilution of 1:1000



Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:107636(TRAF3 Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)