

# transcription termination factor-Specific antibody



Catalog Number: 107637

## Product name

Monoclonal

transcription termination factor-Specific antibody

## Ig Type

## Specificity

Mouse IgG2a

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

transcription termination factor-Specific Mouse Monoclonal antibody. Positive WB detected in HeLa cell. Observed molecular weight by Western-blot: 106 kDa

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_007344). Purification method: Protein A purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107637(TTF1-Specific antibody) at dilution of 1:2000

## Clonality