

# TRIM24 antibody



Catalog Number: 107643

## Product name

TRIM24 antibody

## Immunogen

[Human TRIM24 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TRIM24 Mouse Monoclonal antibody. Positive WB detected in NIH/3T3 cells, HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 116 kDa

## Preparation

This antibody was obtained by immunization of TRIM24 recombinant protein (Accession Number: NM\_015905). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107643(TRIM24 Antibody) at dilution of 1:1000