

# TRIM44 antibody



Catalog Number: 107644

## Product name

TRIM44 antibody

## Specificity

Human; other species not tested.

## Antibody description

TRIM44 Mouse Monoclonal antibody. Positive WB detected in MCF-7 cells, human testis tissue. Positive IF detected in MCF-7 cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of TRIM44 recombinant protein (Accession Number: NM\_017583). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IF: 1:50-1:500

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107644(TRIM44 Antibody) at dilution of 1:1000