

# TSPAN3 antibody



Catalog Number: 107647

## Product name

TSPAN3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TSPAN3 Mouse Monoclonal antibody. Positive WB detected in Hela cells, human blood tissue. Positive FC detected in HeLa cells. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of TSPAN3 recombinant protein (Accession Number: NM\_005724). Purification method: Caprylic Acid/Ammonium Sulfate Precipitation.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107647(TSPAN3 antibody) at dilution of 1:1000



1X10<sup>6</sup> HeLa cells were stained with 0.2ug TSPAN3 antibody (Catalog No:107647, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1000.