

IL-12B antibody

Catalog Number: 107656

Product name

IL-12B antibody

Specificity

Human; other species not tested.

Antibody description

IL-12B Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of IL-12B recombinant protein (Accession Number: NM_002187). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:107656(IL12B antibody) at dilution of 1:5000



-	+	-	+	Overexpression of IL12B-his-myc
+	+	-	-	Rabbit Anti human IL12B polyclonal antibody
-	-	+	+	Mouse Anti his-tag monoclonal antibody

Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107656(IL12B Antibody) at dilution of 1:500