

5 Lipoxygenase antibody



Catalog Number: 107662

Product name

5 Lipoxygenase antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

5 Lipoxygenase Rabbit Polyclonal antibody.
Positive WB detected in HL-60 cells, DU 145 cells.
Positive IF detected in Hela cells. Observed
molecular weight by Western-blot: ~70 kDa

Preparation

This antibody was obtained by immunization of
Recombinant Protein (Accession Number:
NM_000698). Purification method: Protein A
purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH
7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



Immunofluorescent analysis of Hela cells, using
ALOX5 antibody Catalog No:107662 at 1:25
dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

HL-60 cells were subjected to SDS PAGE followed
by western blot with Catalog No:107662(ALOX5
Antibody) at dilution of 1:600