

A2BP1 antibody



Catalog Number: 107665

Product name

A2BP1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

A2BP1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive IHC detected in mouse brain tissue, human brain tissue. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of A2BP1 recombinant protein (Accession Number: NM_145893). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IP


Dilutions

Recommended Dilution:


IP: 1:200-1:2000

IHC: 1:20-1:200

Validations



IP Result of anti-A2BP1 (IP:Catalog No:107665, 4ug; Detection:Catalog No:107665 1:500) with mouse brain tissue lysate 4000ug.



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using Catalog No:107665(A2BP1 Antibody) at dilution of 1:200 (under 10x lens).