

A2LD1 antibody



Catalog Number: 107666

Product name

Polyclonal

A2LD1 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

A2LD1 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse bladder tissue, mouse brain tissue. Observed molecular weight by Western-blot: 17kd

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of A2LD1 recombinant protein (Accession Number: XM_047430708). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:107666(A2LD1 antibody) at dilution of 1:500