

# A4GALT antibody



Catalog Number: 107667

## Product name

A4GALT antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

A4GALT Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells, K-562 cells, Raji cells. Observed molecular weight by Western-blot: 50-52 kDa

## Preparation

This antibody was obtained by immunization of A4GALT recombinant protein (Accession Number: XM\_005261646). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107667(A4GALT Antibody) at dilution of 1:1000