

AARSD1 antibody



Catalog Number: 107676

Product name

Polyclonal

AARSD1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

AARSD1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 45kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of AARSD1 recombinant protein (Accession Number: BC004172). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107676(AARSD1 antibody) at dilution of 1:500