

AATF antibody

Catalog Number: 107678

Product name

AATF antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AATF Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 63kd

Preparation

This antibody was obtained by immunization of AATF recombinant protein (Accession Number: NM_012138). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107678(AATF antibody) at dilution of 1:800



Immunofluorescent analysis of HeLa cells, using AATF antibody Catalog No:107678 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).