

ABCA2 antibody

Catalog Number: 107680

Product name

ABCA2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ABCA2 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, HepG2 cells, L02 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 270 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001606). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:107680(ABCA2 antibody) at dilution of 1:500



Immunofluorescent analysis of (-20oc Ethanol) fixed HepG2 cells using Catalog No:107680(ABCA2 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)