

ABCA8 antibody



Catalog Number: 107682

Product name

Polyclonal

ABCA8 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

ABCA8 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 160-180 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of ABCA8 recombinant protein (Accession Number: bc130280). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107682(ABCA8 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality