

# ABCA8 antibody

Catalog Number: 107682

---

## Product name

ABCA8 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ABCA8 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 160-180 kDa

## Preparation

This antibody was obtained by immunization of ABCA8 recombinant protein (Accession Number: bc130280). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107682(ABCA8 Antibody) at dilution of 1:300