

# ABCB10 antibody



Catalog Number: 107683

## Product name

ABCB10 antibody

## Immunogen

[Human ABCB10 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ABCB10 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells, Jurkat cells. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of ABCB10 recombinant protein (Accession Number: NM\_012089). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107683(ABCB10 antibody) at dilution of 1:500



IP Result of anti-ABCB10 (IP:Catalog No:107683, 4ug; Detection:Catalog No:107683 1:300) with K-562 cells lysate 3200ug.