

# ABCB7 antibody

Catalog Number: 107686

## Product name

ABCB7 antibody

## Immunogen

[Human ABCB7 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ABCB7 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, mouse thymus tissue. Observed molecular weight by Western blot: 83kd

## Preparation

This antibody was obtained by immunization of ABCB7 recombinant protein (Accession Number: NM\_004299). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107686(ABCB7 antibody) at dilution of 1:1000



IP Result of anti-ABCB7 (IP:Catalog No:107686, 4ug; Detection:Catalog No:107686 1:500) with mouse liver tissue lysate 4000ug.