

ACP1 antibody

Catalog Number: 107688

Product name

ACP1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ACP1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue. Positive WB detected in mouse brain tissue, 4T1 cells. Observed molecular weight by Western-blot: 18 kDa

Preparation

This antibody was obtained by immunization of ACP1 recombinant protein (Accession Number: NM_007099). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107688(ACP1 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:107688(ACP1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog

ACP1 antibody



Catalog Number: 107688

No:107688(ACP1 Antibody) at dilution of 1:50 (under 40x lens)