

ACP2 antibody



Catalog Number: 107689

Product name

Polyclonal

ACP2 antibody

Ig Type

Specificity

Rabbit IgG

Human,Mouse,Rat; other species not tested.

Applications

Antibody description

ELISA, WB

ACP2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HepG2 cells, HL-60 cells, mouse brain tissue. Observed molecular weight by Western-blot: 63 kDa,52 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of ACP2 recombinant protein (Accession Number: NM_001610). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107689(ACP2 antibody) at dilution of 1:500