

ACPT antibody



Catalog Number: 107692

Product name

ACPT antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ACPT Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, mouse heart tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of ACPT recombinant protein (Accession Number: NM_033068). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107692(ACPT Antibody) at dilution of 1:1000