

ACPT antibody



Catalog Number: 107692

Product name

Polyclonal

ACPT antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse; other species not tested.

Applications

Antibody description

ELISA, WB

ACPT Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, mouse heart tissue. Observed molecular weight by Western-blot: 35 kDa

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of ACPT recombinant protein (Accession Number: NM_033068). Purification method: Antigen affinity purified.

Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107692(ACPT Antibody) at dilution of 1:1000