

ACRV1 antibody



Catalog Number: 107693

Product name

ACRV1 antibody

Specificity

Human; other species not tested.

Antibody description

ACRV1 Rabbit Polyclonal antibody. Positive WB detected in DU 145 cells, LNCaP cells, PC-3 cells. Positive IF detected in PC-3 cells. Observed molecular weight by Western-blot: 13-36 kDa

Preparation

This antibody was obtained by immunization of ACRV1 recombinant protein (Accession Number: NM_001612). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:20-1:200

Validations



DU 145 cells were subjected to SDS PAGE followed by western blot with Catalog No:107693(ACRV1 antibody) at dilution of 1:4000



Immunofluorescent analysis of PC-3 cells using Catalog No:107693(ACRV1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)