

ACSL3 antibody



Catalog Number: 107696

Product name

ACSL3 antibody

Specificity

human, mouse; other species not tested.

Antibody description

ACSL3 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in mouse kidney tissue, HEK-293 cells. Observed molecular weight by Western-blot: 70-80kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004457).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107696(ACSL3 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using ACSL3 antibody Catalog No:107696 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).