

ACSL5 antibody

Catalog Number: 107698

Product name

ACSL5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACSL5 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, mouse liver tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 66-70kd

Preparation

This antibody was obtained by immunization of ACSL5 recombinant protein (Accession Number: NM_016234). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107698(ACSL5 antibody) at dilution of 1:1200



IP Result of anti-ACSL5 (IP:Catalog No:107698, 4ug; Detection:Catalog No:107698 1:600) with mouse brain tissue lysate 4000ug.