

ACSS1 antibody

Catalog Number: 107702

Product name

ACSS1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACSS1 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in mouse kidney tissue. Observed molecular weight by Western-blot: 75 kDa

Preparation

This antibody was obtained by immunization of ACSS1 recombinant protein (Accession Number: NM_032501). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107702(ACSS1 antibody) at dilution of 1:3000



IP Result of anti-ACSS1 (IP:Catalog No:107702, 4ug; Detection:Catalog No:107702 1:700) with mouse kidney tissue lysate 4000ug.