

ACTL7A antibody



Catalog Number: 107708

Product name

ACTL7A antibody

Immunogen

[Human ACTL7A Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

ACTL7A Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of ACTL7A recombinant protein (Accession Number: NM_006687). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107708(ACTL7A antibody) at dilution of 1:800



Immunofluorescent analysis of HepG2 cells, using ACTL7A antibody Catalog No:107708 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).