

ACBD5 antibody



Catalog Number: 107731

Product name

ACBD5 antibody

Immunogen

[Human ACBD5 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

ACBD5 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 60-70 kDa

Preparation

This antibody was obtained by immunization of ACBD5 recombinant protein (Accession Number: NM_001352579). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107731(ACBD5 antibody) at dilution of 1:500