

ACE antibody



Catalog Number: 107735

Product name

ELISA, WB, IP

ACE antibody

Dilutions

Specificity

Recommended Dilution:

Human, Mouse; other species not tested.

WB: 1:500-1:5000

Antibody description


IP: 1:500-1:5000

ACE Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 170 kDa

Validations

Preparation

This antibody was obtained by immunization of ACE recombinant protein (Accession Number: NM_152831). Purification method: Antigen affinity purified.



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:107735(ACE Antibody) at dilution of 1:1000

Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



IP Result of anti-ACE (IP:Catalog No:107735, 4ug; Detection:Catalog No:107735 1:1000) with mouse heart tissue lysate 4000ug.