

# ACE2 antibody



Catalog Number: 107736

## Product name

ACE2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ACE2 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in mouse kidney tissue, human skeletal muscle tissue. Observed molecular weight by Western-blot: 90-110 kDa

## Preparation

This antibody was obtained by immunization of ACE2 recombinant protein (Accession Number: NM\_021804). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:500-1:5000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107736(ACE2 antibody) at dilution of 1:1000



IP Result of anti-ACE2 (IP:Catalog No:107736, 5ug; Detection:Catalog No:107736 1:1000) with mouse kidney tissue lysate 4000ug.