

ACE2 antibody

Catalog Number: 107736

Product name

ACE2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ACE2 Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Positive WB detected in mouse kidney tissue, human skeletal muscle tissue. Observed molecular weight by Western-blot: 90-110 kDa

Preparation

This antibody was obtained by immunization of ACE2 recombinant protein (Accession Number: NM_021804). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:500-1:5000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:107736(ACE2 antibody) at dilution of 1:1000



IP Result of anti-ACE2 (IP:Catalog No:107736, 5ug; Detection:Catalog No:107736 1:1000) with mouse kidney tissue lysate 4000ug.