

ACOX2 antibody

Catalog Number: 107750

Product name

ACOX2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ACOX2 Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive IF detected in HepG2.

Preparation

This antibody was obtained by immunization of ACOX2 recombinant protein (Accession Number: XM_005265505). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemical of paraffin-embedded human liver using Catalog No:107750(ACOX2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:107750(ACOX2 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using ACOX2 antibody Catalog No:107750 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).