

# ACTR3B antibody



Catalog Number: 107751

## Product name

ACTR3B antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ACTR3B Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HepG2 cells, mouse brain tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 48kd

## Preparation

This antibody was obtained by immunization of ACTR3B recombinant protein (Accession Number: NM\_020445). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107751(ACTR3B antibody) at dilution of 1:800