

ACTR3B antibody



Catalog Number: 107751

Product name

Polyclonal

ACTR3B antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

ACTR3B Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HepG2 cells, mouse brain tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 48kd

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of ACTR3B recombinant protein (Accession Number: NM_020445). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107751(ACTR3B antibody) at dilution of 1:800