

# ACTRT2 antibody



Catalog Number: 107752

## Product name

ACTRT2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ACTRT2 Rabbit Polyclonal antibody. Positive WB detected in rat testis tissue, mouse testis tissue. Observed molecular weight by Western-blot: 40-42 kDa

## Preparation

This antibody was obtained by immunization of ACTRT2 recombinant protein (Accession Number: NM\_080431). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



rat testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107752(ACTRT2 Antibody) at dilution of 1:600