

# ACY3 antibody



Catalog Number: 107756

## Product name

Polyclonal

ACY3 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

ACY3 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 35kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of ACY3 recombinant protein (Accession Number: NM\_080658). Purification method: Antigen affinity purified.

## Validations

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107756(ACY3 antibody) at dilution of 1:400