

ACYP2 antibody



Catalog Number: 107758

Product name

ACYP2 antibody

Immunogen

[Human ACYP2 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

ACYP2 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue. Positive IHC detected in human medulloblastoma. Observed molecular weight by Western-blot: 11kd

Preparation

This antibody was obtained by immunization of ACYP2 recombinant protein (Accession Number: NM_138448). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:107758(ACYP2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human medulloblastoma using Catalog No:107758(ACYP2 antibody) at dilution of 1:50 (under 10x lens)