

# ADAM2 antibody



Catalog Number: 107765

## Product name

Polyclonal

ADAM2 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

ADAM2 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human testis tissue. Observed molecular weight by Western-blot: 50-55 kDa, 65-75 kDa

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of ADAM2 recombinant protein (Accession Number: NM\_001278113). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107765(ADAM2 antibody) at dilution of 1:1000